


## Intelligent Quantifier™ (IQ) for Microsoft® Windows®



*Intelligent Quantifier for Windows is a software system for analyzing scanned images of 1-D and 2-D electrophoresis separations, blots, and colonies/plaques. It allows you to analyze gel files obtained from protein, DNA, RNA, or blot samples and report the data.*

### The Software

- Analyzes protein, DNA, RNA, carbohydrate, colony & blot samples
- Analyzes images from wet gels, dry gels, autoradiograms, blots and photographs
- Allows import from FUJI, Molecular Dynamics and other TIFF-formatted images, and JPEG (JPG) and BMP images

- Displays intensity in optical density (OD), FUJI photo stimulated luminescence (PSL), or Molecular Dynamics units of luminescence
- Allows user to name an individual band, spot, cell, or colony
- Includes an illustrator tool for adding text, arrows, boxes, circles, header and footer information fields onto an image
- Generates band, spot, colony and blot data analysis reports
- Allows user to export images and/or analysis in standard TIFF, Postscript®, and ASCII formats
- Allows user to define the analysis annotation colors and composition
- Includes comprehensive User's Manual hard copy and on-line hypertext
- Provides interoperability with other popular Windows packages (such as Microsoft Word® and Excel®)

### 1-Dimensional electrophoresis analysis

- Automatically detects lanes, finds bands, including irregularly shaped bands, i.e., biphasic bands, "smiling bands," or saturated bands that are substantially wider than the rest of the lane
- Allows the user to vary band detection sensitivity to optimize performance to gel characteristics
- Automatically quantifies bands by defining actual band boundaries based on inflection point analysis
- Subtracts local background

- Provides side-by-side lane profile and band size/quantification comparisons
- Quantifies by whole band and densitometry methods
- Users specify methods of interpolation for standards from three methods: Logarithmic, Linear, and Reciprocal
- Displays single lane density profiles or multiple overlay profiles

### 2-dimensional electrophoresis analysis

- Analyzes very complex (10,000 spots) gels
- Automatically identifies spots in a user-defined area
- Automatically quantifies spots using 36 radians per spot for accurate and reproducible quantification data
- Determines the local background for each spot
- Users specify Logarithmic or Linear methods of interpolation for MW and pI standards

### Colony counting

- Allows colony counting of spots in a user-defined area
- Provides automatic colony identification in a user-defined area
- Analyzes very complex (10,000 colonies) samples
- Determines local background for each colony

### Dot/slot blot analysis

- Flexible enough to accommodate Dot and Slot Blot samples, 1 row by 1 column up to 52 rows by 52 columns
- Provides automatic quantification of individual blot cells
- Includes user-definable background calculation (automatic or manual)
- User can normalize all cells to a specified row, column, or cell
- Concentration standards can be assigned to user-selected cells

## System Requirements

### IBM and compatible systems:

- Processor: Intel® Pentium® or AMD® Equivalent
- Hard Disk Drive: 100 GB (if storing images)
- Operating System: Windows 7, Windows 2000 and Windows XP
- Memory: 2 GB minimum
- USB ports: 2
- Display Resolution: 1280 horizontal on a 14-inch screen

## Ordering Information

### Microsoft Windows for PC

Intelligent Quantifier Full Suite  
(1-D, 2-D, Dot Blot and Colony Count)  
Catalog #100100

### Available Options include:

Full image and analysis systems including compatible image acquisition devices, CPU, Monitor and printer.

To order, contact:

**Bio Image Systems, Inc.**  
721 Bloomfield, Jackson, MI 49203 USA

**Phone: +1.734.604.2623**  
**Fax: +1.517.788.8869**

**email: [info@bioimage.net](mailto:info@bioimage.net)**  
**<http://www.bioimage.net>**

Intelligent Quantifier is a trademark of Bio Image Systems, Inc., Windows, Microsoft, Microsoft Word, and Microsoft Excel are registered trademarks of Microsoft, Inc., Postscript is a registered trademark of Adobe, Inc.