



ACS 5000 and ACS 6000 Specifications

Light Source:

Wavelength	400-750 nm
Color	User-selected Gray, Red, Green, or Blue
Mode	Transmissive and Reflective
Absorbance Range	0.05-3.77 (transmissive), 0-2.0 RD (reflective)
Bit Depth	16 bit (65536 gray levels)
Power	100-240V, 50-60 Hz

Scan Area:

Spatial Resolution	User-selected 100, 150, 300, or 600 dpi
--------------------	---

ACS 5000 (A4)

Sample Area (cm)	26 x 19 (transmissive) 30 x 20 (reflective)
Dimensions (cm)	30 W x 16 H x 50 D

ACS 6000 (A3)

Sample Area (cm)	26 x 19 (transmissive) 30 x 20 (reflective)
Dimensions (cm)	30 W x 16 H x 50 D