

FOR IMMEDIATE RELEASE – July 7, 2003

Contact: Bob Luton, President & CTO
Bio Image Systems, Inc.
721 Bloomfield
Jackson, MI 49203 USA

1-D Electrophoresis Analysis Software Provides 21CFR Part 11 Compliance

Bio Image Systems, Inc. announces the release of *Intelligent Quantifier-11 for Windows 2000*. This software system for analyzing scanned images of 1-D electrophoresis separations is fully compliant with 21CFR Part 11. The software includes fully integrated logins and passwords; record signings (acquire image, save changes, and reports); Analyst, Supervisor, and Administrator access levels; secure audit trails of all operations and a system audit trail.

The *Intelligent Quantifier-11* software analyzes images from wet gels, dry gels, and membranes. It has a complete one-dimensional analysis package including lane definition, band detection, Mw/ IEF standards, and dilution series (concentration standards). The software generates band and lane data analysis reports and comes with an IQ/OQ. Customized reports and complete hardware and software systems are also available.

About Bio Image Systems, Inc.

Bio Image Systems, Inc. is a sales and support company for Bio Image software, including 21 CFR Part-11 Image Analysis Software. Bob Luton, President and Chief Technology Officer, has 18 years experience developing and supporting the Bio Image suite of products. Bio Image products are available for Sun, PC, and Macintosh.

For more information, contact Bob Luton at Bio Image Systems, Inc.

Phone: +1.734.604.2623

Fax: +1.517.788.8869

email: info@bioimage.net

Internet: <http://www.bioimage.net>